## **Analysis appendix**

## **Data Tables**

8cm CCW		8cm CW		
volts (V)	current (A)	volts (V)	current (A)	
148.5	1.355	148.3	1.248	
168.6	1.446	168.4	1.325	
187.9	1.520	188.1	1.400	
207.7	1.604	208.0	1.474	
227.9	1.674	228.2	1.546	
247.5	1.752	247.8	1.613	
10cm CCW 10cm CW		10cm CW		
volts (V)	current (A)	volts (V)	current (A)	
149.0	1.117	148.6	0.960	
168.6	1.177	168.9	1.021	
188.0	1.234	187.6	1.089	
208.2	1.301	208.3	1.165	
227.6	1.367	227.4	1.205	
248.1	1.422	248.0	1.260	

Table 1: Current through Helmholtz coils required to maintain 8 and 10cm beam diameter, respectively, for given anode-cathode voltages (resulting in different beam velocities), in clockwise (CW) and counter-clockwise (CCW) orientations.

## **Variables and Equations**

Variable	e Name	Value and/or SI Unit
$\overline{V}$	Voltage applied to anode	$kg \cdot m^2 \cdot s^{-3}$
I	Current through Helmholtz coils	$C \cdot s^{-1}$
K	Coefficient determining field strength at coil center	$(7.73 \pm 0.04) \cdot 10^{-4} T \cdot A^{-1}$
R	Radius of Helmholtz coils	$0.154 \pm 0.005  m$
$K_R$	Effective flux density for given distance from center of coils	$\delta \alpha T \cdot A^{-1}, \alpha \in [0,1]$
$B_{\scriptscriptstyle T}$	Total effective field in Helmholtz coils	T
$B_{\scriptscriptstyle E}$	Field contributed by Earth's magnetic field	T
$B_l$	Field contributed by current in coil	T
$\boldsymbol{I}_l$	Current in high-current orientation for measurement pair	$C \cdot s^{-1}$
$I_s$	Current in low-current orientation for measurement pair	$C \cdot s^{-1}$
r	Radius of electron beam	m

Table 2: Variables used in all equations and analysis

An electron traveling through a potential difference of V Volts gains kinetic energy

$$eV = \frac{1}{2}mv^2$$

If the electron then travels through a uniform magnetic field  $\vec{B}$ , a force  $e\vec{v} \times \vec{B}$  acts on the electron. Assuming the direction of travel is exactly perpendicular to the field lines of  $\vec{B}$ , we have

$$evB = \frac{mv^2}{r}$$

Rearranging and combining, we get

$$\frac{e}{m} = \frac{2V}{B^2 r^2}$$
.

A series of anode voltage and Helmholtz current data pairs are taken by increasing V, and maintaining the radius r constant. B is calculated using  $K_rI$ . These measurements are repeated for the same voltages with the apparatus rotated  $180^\circ$ , such that the effects of the Earth's magnetic field are counteracted and can be measured.

## Sample calculation

For V = 250, I = 1.6825 (average of CCW and CW measurement pair):

$$\frac{e}{m} = \frac{2.250}{7.73 \cdot 10^{-4} \cdot 1.6825 \cdot 0.08} = 1.847 \cdot 10^{11}$$

The real value is  $1.76 \cdot 10^{11}$ , so this measurement is within 5% of the real value.